



# MED-810A PRO

Q Switch ND YAG Laser Machine



# SUPER PROFESSIONAL

Mini PICO LAEER 5 in 1 laser machine



## Application

Tattoo Removal; PMU Removal;  
Speckle Removal;  
Nevus of Ota Removal;  
Birthmark Removal;  
Skin Rejuvenation, Whiten, Acne  
Tighten Pores



## Application



### 1064nm Tip (Standard)

Professional for dark colors tattoo, pigment removal, nail fungus.



### 532nm Tip (Standard)

Specialized in colorized tattoo treatment, such as red, brown and coffee etc.



### 1320nm Tip (Standard)

Mainly applied to deep cleaning, skin whitening, pore shrinking.



### Focus Tip (Standard) Spot Size 8\*8mm

Specialized for skin whiten & rejuvenation, acne therapy, reduce fine lines, with smooth feeling.



### Fractional Tip (Standard) Spot size 7\*7mm

Specialized for pigmentation issues such as freckles, sunspots, without bleeding and minimal downtime




## Advantages



### 1. TUV Medical CE approved

Strict quality testing process, safe treatment and stable working.



### 2. Aiming light accurate operating

Aiming light can be turn on/off, longer lifetime



### 3. High-efficiency Dermal Pigment & Tattoo Removal

Double energy handpiece especially focus on dermal pigment and colorful tattoo removal, permanent results and less sessions is needed.



### 4. Super long lifetime

Lifetime of xenon lamp more than 1,000,000 shots.



### 5. Intelligent interface

Five different tips for various treatments can be selected

## ■ Technique Specification

---

<b>Model NO.</b>	<b>Med-810A PRO</b>
<b>wavelength</b>	<b>1064nm, 532nm, 1320nm, Fractional, Focus Array; Five Tips Standard</b>
<b>Laser Power</b>	<b>500W</b>
<b>Energy</b>	<b>300-1600mj</b>
<b>Yag bar size</b>	<b>Φ7+Φ8</b>
<b>Display</b>	<b>10.4" True Color LCD Touch Screen</b>
<b>Cooling</b>	<b>Air+Water+Twin Radiators Cooling System</b>
<b>Net Weight</b>	<b>25KG</b>
<b>Dimensions (WxLxH)</b>	<b>40*36*46cm</b>

---