Technique Creates Beauty

MED-340Pro

 $\mathbf{\Sigma}$

ഗ

刀

n


m

Cryolipolysis Fat Reduction Slimming Machine










Application



Big Handle (optional) provides solution for particularly obese parts, removing the fat fastly.



MIDDLE HANDLE

Middle Handle (Standard configuration) gives treatment for whole body, perfect solution for belly, butt and legs.



<

m

 \triangleright

Z

m

S

-

<

⋗

m

SMALL HANDLE 130*60*55 mm


Small Handle (optional) especially offers effective treatment for small area of the body



Technique Creates Beauty

Proven Clinical Results

The equipment adopts advanced cryolipolysis technology. Triglyceride of fat will be converted to the solid (crystalline) at a low temperature of 4-5 degree; Solid triglyceride aging in advance; Normal metabolism makes aging solid triglyceride to be excreted out of body, to get a result of weight loss.



Before and After



After nine weeks of treatment



Before treatment



After 90 days of treatment





Before treatment

After 60 days of treatment

 $\mathbf{\Sigma}$

ഗ

ഗ

m

刀

<

0

m

⋗

Z

σ

m

S


<

⊳

ET I



Before treatment



After 60 days of treatment



Before treatment




After 60 days of treatment


Technique Creates Beauty

Advantages of Cryolipolysis














 Four handles work together Guarantee treatment results and decrease 50% treatment time.

2. Unique Silicone-Contact-cooling Silicone-Contact-Cooling plus Real-Time temperature calibration. Double safe guarantee and more effective. Say goodbye to Metal-Contact-Cooling which with high burning risk.

 Safe Wiring Design
Identified wiring and Plug & Play connector for easy recognition and fast maintenance

4. Intelligent Software, convenient operation Select fat thickness,parameters auto set Four handpieces controlled separately

Independent Modular Design
Modular design ensure easy and convenient
maintenance. Independent water and electrical
circulation avoid short circuit caused by water leakage

6. Double Filtration ProtectionFirst Filtration Protection on back of machine. SecondFiltration Protection on each side of machine. Extendmachine lifetime and ensure convenient maintenance

Technique Specification

Model NO.	Med-340Pro	
Power	1200W	
Focused Technology	Cryolipolysis + Vacuum + LED	
	Totol 4;	
Handpiece	2 x Middle (170*60*55mm)	
	1 x Big (210*60*55mm)	
	1 x Small (130*60*55mm)	
Cooling Stand-by Working	-10°C~0°C Continuous contact cooling	
Vacuum pressure	Level 1~10 adjustable	
Display	10.4" True Color LCD Touch Screen	
Net Weight	68 KG	
Dimensions	480*420*1050mm	

⊳

Z

m

S

-1

<

⊳

-

m